FR

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EN

**CIMPAT L.C.** is a light-curing material for temporary filling and sealing. The elastic consistency assures an easy removal, even in case of undercuts, and finishing the cavity is not necessary. Thus the fitting of the temporary filling will be improved.

#### Indications

All types of temporary fillings, temporary protection in inlay and onlay technique.

#### Application

Preparation, lining and impression to be carried out as usual. Rinse and dry the cavity.

Apply **CIMPAT L.C.** into the cavity with a suitable instrument (round condenser, spatula, etc.).

#### Further handling

Contour with slight excess. For use in the inlay and onlay technique, keep the approximate contact points free of resin with dental floss or a matrix before curing to maintain tooth mobility.

Cure for 20 - 40 s (depending on layer thickness) with a conventional polymerisation device. The light output should be a minimum of 500 mW/cm<sup>2</sup> on halogen polymerisation devices and 300 mW/cm<sup>2</sup> on LED devices. As a rule, 5 mm deep layers are sufficiently cured in 40 s. Increased exposure times are required when the light source cannot be brought close enough to the material.

Check occlusion. Excess can be removed with a scalpel. Silicone points are suitable for grinding.

#### Removal of the temporary filling

Insert a rigid instrument into the temporary filling and extract it.

#### Advice and precautions

Adhesion to dental resins or linings (e.g. light-curing materials) may occur. In such cases apply a suitable insulating agent (e.g. glycerine, vaseline). Avoid use of lining materials containing eugenol, since these could impede the setting of **CIMPAT L.C.**

In case of more than two weeks' stay in the mouth, and in occlusal or marginal areas, reinforcement with composite is recommended. Apply a small quantity of a composite material onto the corresponding areas and polymerise after contouring together with **CIMPAT L.C.**

The slight expansion of **CIMPAT L.C.** ensures a tight seal of the temporary restoration. When treating adjacent teeth simultaneously, remove excess material in approximal areas prior to polymerisation to maintain tooth mobility and contact points.

**CIMPAT L.C.** contains hydroxyethylmethacrylate, BHT, acrylate ester and polymers.

Allergies to ingredients are possible. Avoid contact of unpolymerised paste with mucous membranes to avoid sensitivity reactions. Do not use after expiry date.

#### Storage

Store at temperatures of 4°C to 23°C. Refrigeration is advised if the material is not used regularly.

Our preparations have been developed for use in dentistry. As far as the application of the products delivered by us is concerned, our verbal and/or written information has been given to the best of our knowledge

and without obligation. Our information and/or advice do not relieve you from examining the materials delivered by us as to their suitability for the intended purposes of application. As the application of our preparations is beyond our control, the user is fully responsible for the application. Of course, we guarantee the quality of our preparations in accordance with the existing standards and corresponding to the conditions as stipulated in our general terms of sale and delivery.

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