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#### PACKAGING AND STERILISATION

- Bone substitutes are delivered sterile.
- They are sterilised with gamma radiation at a minimum dose of 25 kGy.
- The integrity of the individual sterility protector should be checked before use.
- Do not use if the packaging is damaged ☷

#### RE-STERILISATION AND REUSE

The re-sterilisation and/or reuse of RTR bone substitutes is strictly forbidden, as this may expose the patient to contamination that may lead to infection. The re-sterilisation and/or reuse of a RTR bone substitute would result in a decrease in the performance of the device, namely reduced mechanical strength (risk of fracture) and an incomplete bone remodelling phenomenon. There is also a risk of a significant inflammatory reaction (for example, a host reaction against the graft if it has been implanted before).

#### ELIMINATION AFTER USE

In the case of extraction of the device, please handle the components in accordance with the centre's medical waste disposal procedure.

#### INDICATIONS

RTR bone substitutes are intended for filling and reconstructing bone loss.

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Before using a RTR bone substitute, the surgeon must take cognisance of the following:

#### GENERAL INFORMATION

The manufacturer of RTR bone substitutes guarantees the quality of manufacture and material.

Several criteria are of decisive importance for the success of an intervention using a bone substitute:

- The appropriate selection of the patient for the intervention in question
- The strictest asepsis

The surgeon is responsible for the complications that may result from an erroneous indication and a faulty operating technique or fault involving asepsis, and these cannot be imputed to the manufacturer.

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#### MANUFACTURING - MATERIALS

The RTR bone substitutes are made of tricalcium phosphate  $\beta$  (TCP;  $\text{Ca}_3(\text{PO}_4)_2$ ; origin: synthetic) according to the latest technologies to obtain an excellent level of purity. However, their mechanical and functional characteristics can only be guaranteed if the instructions in this manual have been observed.

#### CONTRAINDICATIONS

The following circumstances are a contraindication to the placement of grafts:

- General contraindications against grafts
- Acute osteomyelitis or necrosis of the recipient site
- Bone degeneration
- Known active infection
- Allergy to a constituent of the bone substitute
- When any surgery is contraindicated

#### STORAGE AND HANDLING

- Bone substitutes must be stored in their original packaging.

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#### COMPLICATIONS

The following complications have sometimes been observed, and these require special attention:

- Delayed healing
- Fracture
- Delayed union and pseudarthrosis
- Pain
- Temporary paralysis of a nerve
- Deep vein thrombosis and pulmonary embolism
- Infection
- Inflammation
- Oedema of the wound
- Allergic reaction
- Opening of the wound
- Reaction to a foreign body
- Tissue degeneration
- Necrosis
- Mobilisation of the dental implant held in place (in whole or in part) by the bone substitute
- Benign nosebleeds (benign epistaxis)
- Benign positional paroxysmal vertigo
- Postoperative seroma
- Failure of bone grafting
- Non-resorption of the bone substitute

#### RE-STERILISATION AND STERILISATION

- Bone substitutes are delivered sterile.
- They are sterilised with gamma radiation at a minimum dose of 25 kGy.

- The integrity of the individual sterility protector should be checked before use.

- Do not use if the packaging is damaged ☷

#### PACKAGING AND STERILISATION

- Bone substitutes are delivered sterile.
- They are sterilised with gamma radiation at a minimum dose of 25 kGy.

- The integrity of the individual sterility protector should be checked before use.

- Do not use if the packaging is damaged ☷

#### RE-STERILISATION AND REUSE

The re-sterilisation and/or reuse of RTR bone substitutes is strictly forbidden, as this may expose the patient to contamination that may lead to infection. The re-sterilisation and/or reuse of a RTR bone substitute would result in a decrease in the performance of the device, namely reduced mechanical strength (risk of fracture) and an incomplete bone remodelling phenomenon. There is also a risk of a significant inflammatory reaction (for example, a host reaction against the graft if it has been implanted before).

#### ELIMINATION AFTER USE

In the case of extraction of the device, please handle the components in accordance with the centre's medical waste disposal procedure.

#### INDICATIONS

RTR bone substitutes are intended for filling and reconstructing bone loss.

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#### PERFORMANCES ATTENDUES

The surgeon must choose the size and shape of the graft best suited to the surgical treatment. The radiological examination will allow him to make this choice.

The bone site intended to be in contact with the graft must be abraded and, if possible, well vascularised in order to have a good quality bone bed.

In order for the graft to function properly, it must be placed in close and stable contact with the bone site (by compression or osteosynthesis system). In addition, to improve the graft, it should be well recommended to impregnate the implant

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