**FR**

## Composite fluide N’Durance® Dimer Flow

Technologie de conversion Nano-Dimer™

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**UK**

## Flowable Composite N’Durance® Dimer Flow

Nano-Dimer conversion technology™

<span></span>
Low-Shrinkage
Nano-Hybrid Flowable Composite
Light-Curing, Low Viscosity
For Professional Dental Use.

**FORMULATION :**

N’Durance Dimer Flow flowable composite is a 41% filled (by volume) composite based on a tri-monomer system which includes dimer-acid dimethacrylate. The filler particle sizes range from a nanosized ytterbium trifluoride to a submicron barium glass. It contains 2% sodium fluoride (by weight), and contains about 1wt% catalysts, inhibitors and pigments.

**INDICATIONS :**

N’Durance Dimer Flow flowable composite is indicated for class III and V restorations. It can also be used for repair of crown margins at the discretion of the dentist. In posterior composite restorative techniques N’Durance Dimer Flow may also be used immediately after the bonding adhesive as a liner at the base of the restoration and to modify the margins of (bisacryl) temporary crowns. N’Durance Dimer Flow provides solutions to a number of clinical problems in dental restorative technique especially where flow can improve adaptation.

**PROPERTIES :**

N’Durance Dimer Flow flowable composite is a light-curing, resin based, dental restorative material. It is available in a selection of shades matched to the VITA® Shade Guide. Its low viscosity and precise syringe delivery make it flowable, easy to place and finish, and adaptable to a variety of clinical situations.

**INSTRUCTIONS FOR USE :**

- The tooth should be cleaned by scaling and prophylaxis with flour of pumice.
- Shade selection: Esthetic choices should be determined prior to isolation and/ or preparation of the tooth/teeth. Do not stare at the shade tabs and teeth for more than 10-15 seconds during shade selection. Use of a VITA Shade Guide is recommended.
- Isolation: Isolate the teeth. Use of a rubber dam is highly recommended.
- Cavity preparation: Follow usual procedures for tooth cavity preparation, leaving no residual material or base from any previous restoration.
- Pulp Protection: Base the preparation as needed. Calcium hydroxide cement may be used for covering small pulpal exposures. Any additional dentinal coverage desired can be achieved by using a glass ionomer lining cement. Eugenol containing cements must be avoided.
- Bonding: Use a state of the art dentin/enamel bonding system to bond N’Durance Dimer Flow composite to the tooth structure. Other modified bonding systems can be used to bond N’Durance Dimer Flow to other kinds of materials used in restorative techniques. Cure the bonding resin before placing the composite. Follow the manufacturer’s instructions.
- Dispensing N’Durance Dimer Flow: Attach a syringe tip to the N’Durance Dimer Flow flowable composite syringe. N’Durance Dimer Flow can be applied directly to the tooth surface by exerting gentle pressure on the syringe plunger. Stop pressure on plunger before the desired amount of material is dispensed, since the flow of material from the syringe tip continues for a short time after the pressure is released. Slight back pressure on the plunger can help to decrease the continued flow of the resin out of the syringe.
- Placement and Curing: A layer of N’Durance Dimer Flow should not exceed 2 to 2 mm. Light cure each layer for 30 seconds with a blue, dental curing light with a wavelength around 470-480 nanometers and a minimum of 400mW/cm intensity. For Universal Opaque shade, apply one thin coat ( mm thick) to cover the discolored or metal area only and cure as previously stated.
- Finishing: After curing, remove flash with a sharp instrument. N’Durance Dimer Flow provides excellent smoothness and may simplify finishing. If needed, finish using a fine diamond, multifluted carbide finishing bur, or finishing disc and strips.
- Polishing: Various polishing materials are available that give excellent results.

**WARNINGS AND PRECAUTIONS FOR USE :**

- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection
- N’Durance Dimer Flow flowable composite contains dimethacrylate resins. Avoid use of this product on patients with known dimethacrylate allergies. To reduce the risk of allergic response, avoid exposure to uncured resin.
- Do not use eugenol-containing materials for pulp protection since they can retard the curing process.
- Do not allow saliva or water contamination of etched tooth surface or unset material.
- Any non-sterilizable item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.

**Xi - Irritant**

- Irritating to eyes and skin
- Contact with acids liberates very toxic gas
- In case of contact with eyes, rinse immediately with plenty of water and seek medical advice
- Wear suitable gloves
- This material and/or its container must be disposed of as hazardous waste
- Refer to special instructions/safety data sheets.

**STORAGE :**

- This composite is designed to be stored between 10°-24°C.
- Replace caps immediately after dispensing.
- If stored in a refrigerator, bring to room temperature prior to use (21°C - 24°C).
- Do not store material in proximity of eugenol-containing materials.
- Do not expose to elevated temperatures or intense light.

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Vita® is a trademark of VITA Zahnfabrik, Bad Säckingen, Germany.