**N'Durance®**

Technologie de conversion Nano-Dimer®

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FORMULATION:
The unique resin based matrix contains approximately 19 wt % of ethoxylated BisGMA, UDMA and the new dicarbamate dimethacrylate dimer acid.

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DIRECTIONS FOR USE:
1. Teeth should be cleaned by scaling and prophylaxis with flour of pumice.

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WARNINGS:

- Avoid contact with skin and eyes.
- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.
- N'DURANCE contains dimethacrylate resins. Avoid use of this product on patients with known dimethacrylate allergies. To reduce the risk of allergic response, avoid exposure to uncured resin.

• Do not use eugenol-containing materials for pulp protection since they can retard the curing process.

• Do not allow saliva or water contamination of etched tooth surface or unset material.

• Any non-sterile item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.

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N'Durance®

Nano-Dimer conversion technology®

Light-curing, low-shrinkage, multi-purpose, esthetic, nanohybrid dental restorative material.

For professional dental use.

INDICATIONS:

N'DURANCE is indicated for direct esthetic restoration of anterior and posterior teeth. It provides exceptional esthetics since it is formulated in shades matching the VITA® Shade guide. This low shrinkage nanohybrid composite can be used in all classes of direct dental restorations, Class I, II, III, IV and V.

DESCRIPTION:

- N'DURANCE is a resin based light-cured dental restorative material with low volumetric contraction. The low shrinkage produces less stress on the bond with the tooth.
- N'DURANCE has a high polymerization conversion. This high degree of conversion is achieved with short curing times and creates improved physical and mechanical properties, and less elution of any residual monomer.
- N'DURANCE is a high viscosity condensable composite that will resist sticking to instruments, yet will have sufficient flow to adapt to the walls of the cavity preparation.
- N'DURANCE contains a yttrium fluoride mono-dispersed nano-sized particle which combined with regular silica and barium glass composes a state of the art nano-filler system.
- N'DURANCE is radio-opaque and provides high compressive and flexural strengths along with a low water uptake.

FORMULATION:

The unique resin based matrix contains approximately 19 wt % of ethoxylated BisGMA, UDMA and the new dicarbamate dimethacrylate dimer acid.

The filler system contains approx. 80 wt% (65 vol%) silanated 40 nanometer ytterbium fluoride, silanated 500 nanometer barium glass and 10 nanometer silica.

There is approximately 1 wt% of catalyst, inhibitors and pigments.

DIRECTIONS FOR USE:

1. Teeth should be cleaned by scaling and prophylaxis with flour of pumice.

2. Shade selection should be performed prior to isolation and/or preparation of the teeth. Do not stare at the tabs and teeth for more than 10-15 seconds during shade selection. Use of a VITA shade guide is recommended.

3. Isolate the teeth. Use of a rubber dam is highly recommended.

4. Follow usual dental procedure for cavity preparation, leaving no residual material or base from any previous restoration.

5. Base the preparation as needed. Calcium hydroxide cement may be used for covering small pulpal exposures. Any additional dentinal coverage desired can be achieved by using a glass ionomer lining cement. Eugenol containing cement must be avoided.

6. Place the appropriate matrix or crown form to assure proper proximal contour of the tooth. Place wedges to produce good gingival adaptation and teeth separation. Burnish the matrix to achieve proximal contact. The matrix may be placed following the adhesive application step if preferred.

7. Use a state of the art dentin/enamel bonding system to bond this composite to the tooth structure. Modified bonding systems can be used to bond this composite to other kinds of materials used in restorative techniques. Cure the bonding resin before placing the composite. Follow the manufacturer's instructions.

8. Dispense the composite following the directions corresponding to the dispensing system chosen:

Syringe: If the product is packaged in a syringe, dispense onto a mixing pad the amount of composite needed to complete the restoration(s). Protect the dispensed material from direct light exposure with an opaque cover. Avoid placing under direct intensity of the operator light during the delivery process.

Single Dose Capsules: Insert the capsule into a dispenser. Deliver an increment of the composite directly into the preparation.

9. The composite should be placed into the preparation in increments no greater than 2.5 mm. It is important to use thinner layers for darker shades. Place the composite into the cavity using the appropriately designed dental hand instrument and adapt the composite to the walls of the preparation. Should the material begin to stick to the instrument, the instrument can be wiped clean with a dry, alcohol gauze.

10. Using any acceptable dental curing lamp, cure each layer for 30 seconds (we recommend using a quartz-tungsten-halogen curing light with a minimum radiance of 500mW/cm² in the spectral range of 470-480 nanometers). An additional 20 seconds may be used from other directions.

11. After curing, remove flash with a sharp instrument. If necessary adjust restoration to proper margins, contours, contacts and occlusion. Finish using a fine diamond, multifaceted carbide finishing bur, or finishing disc and strips.

12. Rubber cups and points and various polishing pastes are available which can enhance the polish.

13. After the polishing procedure is completed, microscopic surface defects can be resealed. The composite surface can be cleaned and the surrounding enamel re-etched with phosphoric acid, washed and dried thoroughly. The bonding agent is applied in a thin layer, dried and cured.

WARNINGS:

- Avoid contact with skin and eyes.
- Dentists and assistants should wear gloves and protective eye wear. Patients also should wear eye protection.

• N'DURANCE contains dimethacrylate resins. Avoid use of this product on patients with known dimethacrylate allergies. To reduce the risk of allergic response, avoid exposure to uncured resin.

- Do not use eugenol-containing materials for pulp protection since they can retard the curing process.

- Do not allow saliva or water contamination of etched tooth surface or unset material.

- Any non-sterile item that is handled in the delivery of the dental service should be disinfected by standard dental office hygiene procedures.

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VITA® is a registered trademark of VITA Zahnfabrik, Bad Säckingen, Germany.