FR

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EN

Dental Injection Needles, with patented scalpel bevel
Sterile, siliconized, single use, with bevel mark

INSTRUCTIONS FOR USE

- Break the sterility guarantee seal just before use by turning the two parts of the protective sleeve in opposite directions.
- Screw the pre-threaded needle onto the pre-threaded hub of the syringe.
- Remove the plastic cap just before the injection is administered. It can be used again to remove the needle from the syringe to avoid risk of needle stick injury.
- Discard the used needle into a puncture-resistant container in accordance with all applicable medical waste regulations.

WARNINGS

- Never recap needles by hand. Use only accepted dental recapping procedures. Recapping needles by hand or recapping bent or broken needles can cause unintended needle sticks.
- Do not bend, break, or otherwise stress needles as serious injuries to you and/or your patient can occur.
- Do not insert needle to hub during injections, as needles can break and become lodged in patient's tissue, potentially causing serious permanent injury. Do not use short needles when the expected depth in soft tissue

approaches the length of the needle.

- Avoid extreme pressure and excessive needle movement during injection as this may cause needles to break, which may result in serious injury to you and/or your patient.
- Be especially vigilant when repositioning a needle in a patient who appears to be apprehensive.

CAUTIONS

- Store in original packaging until ready to use.
- Read the Instructions For Use supplied with your injection device.
- If seal is broken, the needle is NOT STERILE and MUST NOT BE USED.
- For single use only: This device is for use in a single patient. Any reuse exposes other patients and the user to critical risks, including infection and trauma since the needle is no longer sterile and is more likely to break.
- If multiple injections are required, replace the needle for each cartridge in order to minimize the risk of needle point damage.
1 patient / 1 cartridge / 1 needle

PRESENTATION

100 two-point disposable non-sterilizable dental needles

For professional dental use only

AU - Sponsor: Specialities Septodont Pty Ltd
 14 Cliffbrook Crescent, Emu Plains, NSW 2750
 Australia

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